

NEITZ

NEITZ BINOCULAR LOUPE BL β -TITANIUM FRAME SERIES

- Newly developed β -Titanium made frame is Durable and Non-breakable.
- The spring nose support and the forehead pad realize comfortable wearing even in the long time use.
- The adjustment by the flexible loupe hinge accommodates the loupe position in all the angle and distance for acquiring the best view.
- Hinge stopper assists the recovery of the position of loupe, when flipping down the loupe hinge.



BLS-2 New Model Debut !




Serves bright and wide field of view, suitable for cardiovascular surgery

Made in Japan

NEITZ BINOCULAR LOUPE BL β -TITANIUM FRAME SERIES


New Feature

Spring Nose Support (Patent pending)



Nose pad is associated with a spring that softens touch to the nose and hold the loupe frame firmly.

Flip-Up Stopper



Hinge stopper assists to recover the original position of loupe, when flipping it down.

Specification

Model of Loupe	Galilean				Prismatic	
	BLS-1	BLS-2	BLS-3	BLD-3	BLP-4	BLP-6
Magnification	1.5x-2.0x	2.0x-2.5x	2.5x-3.0x	2.5x-3.0x	4.0x-4.5x	5.5x-6.0x
Working distance (mm)	550-350	550-350	550-350	350-250	550-350	550-350
Field of view (mm)	ϕ 145-120	ϕ 113-58	ϕ 78-38	ϕ 80-50	ϕ 80-50	ϕ 40-25
Frame color	4 colors : RED, BLUE, GOLD & BLACK are available with your choice of magnifications.					



Red (with BLP-4)



Blue (with BLS-2)



Gold (with BLP-4)



Black (with BLS-3)

Accessories and Option



LED Spot Illuminator NSI-III

- High-Powered LED provides Distortion-free, Clean-cut Spot Illumination
- Maximum Intensity(14000Lx) : 4 hours (at 40 cm from front of illuminator)
- Intensity and Spot Size Adjustable



BLP-4 with NSI-III, Side Shield and Strap (NSI-III is option)

NEITZ

NEITZ INSTRUMENTS CO., LTD.
 36-13 Wakamatsu-Cho, Shinjuku-Ku, Tokyo 162-0056, Japan
 Phone: (03) 3204-2068 Fax: (03) 3204-2069
 E-Mail : int@neitz.co.jp
 URL : www.neitz-ophthalmic.com